

Analysis Date:

Re-test Date:

ADAGRASIB IMPURITY 1

Identification


Chemical Name : 2-((2S)-4-(7-(2,3-Dimethylphenyl)-2-(((S)-1-methylpyrrolidin-2-yl)methoxy)-5,6,7,8-tetrahydropyrido[3,4-d]pyrimidin-4-yl)-1-((E)-4-(3-fluoropyrrolidin-1-yl)but-2-enoyl)piperazin-2-yl)acetonitrile

CAT No : ALL-A07803

CAS No : 2326521-66-6

Molecular Formula : C₃₅H₄₇FN₈O₂

Molecular Weight : 630.8



Analytical Information

Batch Code : ALL-A07803

HPLC Purity : 98.00 %

Solubility : MeOH: ACN

Potency : N.A

Appearance of Product : Off White Solid

Mass : Confirm

Long Term Storage : 2-8°C

IR Analysis : Confirm

Weight Loss By TGA : N.A

¹H NMR : Confirm

Residue Of Ignition : N.A

Additional Information

%Potency = [100 -(Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 = N.A

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Swapnil Kausatkar Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, ¹H NMR, IR

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E), Mumbai, Maharashtra 421203