

Analysis Date: 23-05-2025

Re-test Date: 23-05-2028

SALBUTAMOL 4-O-SULFATE

Identification


Chemical Name : 4-(2-(Tert-butylamino)-1-hydroxyethyl)-2-(hydroxymethyl)phenyl hydrogen sulfate

CAT No : ALL-A08180

CAS No : 63908-03-2

Molecular Formula : C₁₃H₂₁NO₆S

Molecular Weight : 319.4



Analytical Information

Batch Code : ALL-A08180

HPLC Purity : 98.06%

Solubility : MeOH: ACN

Potency : 96.70%

Appearance of Product : Off White Solid

Mass : Confirm

Long Term Storage : 2-8 °C

IR Analysis : Confirm

Weight Content By TGA : 0.875%

¹H NMR : Confirm

Residue Of Ignition : 0.509%

Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 - (0.875+0.509 × 98.06)]/100 = 96.70%

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse Jr. Executive	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, ¹H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E), Mumbai, Maharashtra 421203