


TENOXICAM D3

Identification

Chemical Name : Tenoxicam D3
CAT No : ALL-D09313
CAS No : 2140327-41-7
Molecular Formula : C₁₃H₈D₃N₃O₄S₂
Molecular Weight : 340.4



Analytical Information

Batch No	: ALL-D09313	HPLC Purity	: 98.10 %
Solubility	: MeOH: ACN	Potency	: 96.25 %
Appearance of Product	: Yellow Solid	% D3	: 99.00 %
Long Term Storage	: 2-8 °C	Mass	: Confirm
Weight Loss By TGA	: 1.323 %	IR Analysis	: Confirm
Residue Of Ignition	: 1.332 %	¹H NMR	: Confirm

Additional Information

%Potency = [100 - (Weight Loss By TGA % + Residue Of Ignition %) × Chromatographic Purity%]/100 =
 [100 - (1.323+1.332) × 98.10]/100 = 95.49 %

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse QA Officer	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, ¹H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature, unless
 otherwise specified

Corporate Office